MR-Guided Intervention & Therapy on the GE MR Platform

Aaron Flammang Product Manager – 1.5T and software Jan. 28, 2016

Interventional CMR Workshop Los Angeles




GE MR Surgical Suite






2-room Suite Installed Base of 30+

Role of MR in CMR Intervention / Therapy




Triage & Pre-Therapy Planning

- MR sequences to assess anatomy, function, flow, lesions
- Patient-specific surgical planning, e.g. mitral valve, TAVI
- 3D Printing from High-Res 3D MR



Procedural Image Guidance

- Real-time catheter tracking
- Therapy Monitoring: Thermometry, T1/T2 Mapping, Late Enh




Post-Therapy Follow-up

• MR sequences to assess therapy & guide medical treatment




New Application Development Relevant to Cardiac MR-Guided Therapy


ViosWorks



4D Flow



Delayed Enhancment




High-Res 3D Late Enhancement




Robust Single Shot MDE

T1 & T2 Mapping



New Application Development Relevant to Cardiac MR-Guided Therapy

External Scan Interface - ExSI



Enables remote real-time control for third parties using standard interfaces



GE Interventional MR Collaborators

Academic Site	Project
Sunnybrook	Cardiac EP: RT guidance, visualization,
University of Dundee	In-room scan control
Stanford	Needle tracking; MR thermometry
UCSF, Zurich Kinderspital	MRgFUS
University of Wisconsin	Peripheral vascular catheter
Robarts	Cardiovascular catheter tracking
Industry	
HeartVista	Real-time scan control
InseRT MR	Software for device control
InSightec	MRgFUS technology















